

# InTouch Vita<sup>®</sup>

## User Guide

# Table of Contents

COPYRIGHT .....	5
SAFETY INSTRUCTIONS .....	7
User Profile .....	7
Notes, Cautions, and Warnings .....	7
Safety Symbols .....	7
Electromagnetic Compatibility - Guidance and Manufacturer's Declaration .....	9
Table 1: Electromagnetic Emissions .....	10
Table 2: Electromagnetic Immunity .....	11
Table 3: Electromagnetic Immunity .....	13
Table 4: Recommended separation distances .....	15
Safety Warnings and Cautions .....	16
VIRTUAL CARE SYSTEM OVERVIEW .....	19
INTOUCH VITA AND ITS INTENDED APPLICATION .....	21
InTouch Vita Indications for Use .....	21
HIPAA POLICIES AND PROCEDURES .....	22
Recommended Guidelines for HIPAA Compliance .....	22
Access to Provider Access Software .....	22
Discussion and Display of PHI .....	23
Stored Images and Video .....	23
Disclosure of PHI .....	24
INTOUCH VITA OVERVIEW .....	25
INTOUCH VITA BASICS .....	28
Charging the InTouch Vita .....	28
Charging Status Indicator Lights .....	28
Vita Power Controls .....	29
Power On Sequence .....	29
Power Off Sequence .....	30
Vita Emergency Stop .....	31
To disable Vita movement: .....	31
To enable Vita movement again: .....	31
Vita Volume Adjustment Knobs .....	32
Vita Head and Chest Display .....	33
Vita Touch-Screen Chest Display .....	34
Vita AutoDRIVE Screen .....	35
Docking the Vita .....	36

Change AutoDRIVE Map .....	36
Vita Settings Screen .....	38
InTouch Vita Wi- Fi Connection Setup .....	39
Connectivity Details .....	40
Device Cleaning Mode .....	42
Date and Time Settings .....	43
Device Settings - More Options .....	44
Vita Change Map .....	45
Vita Edit Badge .....	46
Vita Status Indicator Lights .....	48
CLEANING AND MAINTENANCE .....	50
Cleaning Instructions .....	50
Approved Disinfectants .....	51
Cleaning Sensor Windows (Vita) .....	52
Maintenance and Inspection .....	52
InTouch Vita Maintenance .....	52
Vita Inspecting the ESD Cable .....	52
VITA STETHOSCOPE .....	54
Vita Stethoscope Components .....	54
Vita Stethoscope Operation .....	54
Vita Stethoscope Headset .....	54
Technical Service—Stethoscope .....	54
PRIVACY HANDSET .....	56
INTOUCH VITA AUX VIDEO PORTS .....	57
INTOUCH VITA PHYSICAL SPECIFICATIONS .....	58
Vita Technical Specifications .....	59
Vita System Input Power Requirements .....	60
Vita Classification.....	60
<b>Vita Environmental Specifications .....</b>	<b>60</b>
<b>Operating .....</b>	<b>-60</b>
<b>Non-operating/Storage/Transport .....</b>	<b>60</b>

NETWORK CONFIGURATION .....	61
Configuring the InTouch Patient Access Device Wireless Connection .....	61
Non-Overlapping Channels .....	61
Transmitting Power .....	61
Interference .....	61
Security Options .....	62
Network Installation .....	63
Network Requirements .....	64
Provider Access User Authentication .....	64
Firewall Requirements .....	64
Video Information .....	64
Bandwidth Requirements .....	65
Line Quality Requirements .....	65
Wireless Network Requirements .....	65
Satellite Networks .....	66
Encryption .....	66
Virus Protection .....	66
CONTACT INFORMATION .....	67
24/7 Live Technical Support .....	67
24/7 Live Remote Assistance & 24/7 Live Chat .....	67
E-mail Support .....	67
Website .....	67
Sales & Product Demos .....	67

# Copyright

© 2019 InTouch Technologies, Inc. All rights reserved. InTouch Vita and this manual contain information including, but not limited to, instructions, descriptions, definitions, firmware and software, which are proprietary to InTouch Technologies, Inc. This copyrighted information is protected under Title 17 of the United States Code. This information shall not be copied, modified, or used in any manner which may violate any rights of InTouch Technologies, Inc. InTouch Technologies, Inc. will strictly enforce all of its rights to this information.

Patent(s): <http://www.intouchhealth.com/patents>

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/ or other materials provided with the distribution.
- Neither the name of the Xiph.org Foundation nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

Speex Open Source Software includes:

Copyright 2002-2008	Xiph.org Foundation
Copyright 2002-2008	Jean- Marc Valin
Copyright 2005-2007	Analog Devices Inc.
Copyright 2005-2008	Commonwealth Scientific and Industrial Research Organisation (CSIRO)
Copyright 1993, 2002, 2006	David Rowe
Copyright 2003	EpicGames
Copyright 1992-1994	Jutta Degener, Carsten Bormann

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A

PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE FOUNDATION OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

Microsoft® and Windows® are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Last updated: December, 2019

# Safety Instructions

## User Profile

Trained healthcare professionals are the intended users of any InTouch Health Patient Access Device. Users of the system require clinical judgment and experience to review and interpret the patient data transmitted.

## Notes, Cautions, and Warnings

Information in **ORANGE BOLD FACED TYPE** throughout this User Guide should be interpreted in the following context:

**NOTE:** Provides supplementary information for facilitating operation of the system.

**CAUTION:** Presents instructions for avoiding damage to the system.

**WARNING:** Disregarding this information may prove hazardous to the safety of a person near the InTouch Health Patient Access Device.

## Safety Symbols

Symbols appearing on labels on the Patient Access Device and other equipment include the following:

<b>Attention</b> – Consult accompanying documents for description of intended use.	
<b>Warning Dangerous Voltage</b> – Touching exposed contacts may cause electrical shock. Safety features designed into the device do not allow exposed live AC contacts until the Patient Access Device is fully engaged with its docking station. When fully engaged, contacts are not accessible.	
<b>Wireless Transmitter Notification</b> – Non-ionizing electromagnetic radiation. This device communicates over the 802.11 ac/ a/ b/ g/ n standard for wireless communication.	

<p><b>Body Float</b>—Type B. A Type B piece of equipment is one that provides a particular degree of protection against electric shock, particularly regarding allowable leakage current and reliability of the protective earth connection (grounding).</p>	
<p><b>Body Float</b> – Type BF applied part. Type B equipment with an F-Type applied part. An F-Type applied part is one that extends from the patient into the equipment and is isolated from all other parts of the equipment.</p>	
<p><b>Consult Operator's Manual</b> – Operating Instructions are contained in a separate instruction manual.</p> <p>Some Patient Access Devices utilize a Class II laser which complies with 21 CFR Chapter 1, subchapter j. Maximum laser radiation output is less than one milliWatt. Refer to the appropriate Patient Access Device User Manual to determine the specific wavelength (color) employed by the laser pointer.</p>	



The InTouch Vita utilizes a Class II Laser which complies with 21 CFR Chapter 1, Sub-Chapter J. Maximum Laser radiation output is less than one milliwatt at a wavelength of 635 nanometers.

**CAUTION:** Dazzle, flash-blindness, and afterimages may be caused by a beam from a Class II laser product, particularly under low ambient light conditions. This may have indirect general safety implications resulting from temporary disturbance of vision or from startle reactions. Such visual disturbances could be of particular concern when performing safety-critical operations. Provider Access Software users should not direct the laser beam at persons or at reflective surfaces that may cause disturbances.

## Electromagnetic Compatibility - Guidance and Manufacturer's Declaration

The InTouch Vita system complies with IEC 60601-1-2 , General Requirements for Safety—Collateral standard: Electromagnetic compatibility. Performance of the device is unaffected by exposure to the compliance levels described in Tables 1,2,3 and 4 in the following section.

Special precautions and installation information for the InTouch Vita for electromagnetic compatibility (EMC) are provided below:

- Equipment in hospital environments, including the InTouch Vita and other portable or mobile communications equipment, can produce Electromagnetic Interference (EMI), which may affect the function of these devices. Such effects are prevented by use of equipment with EMI characteristics proven below recognized limits, as identified in the tables below.
- In the event of suspected interference from other equipment, which prevents the proper functioning of the InTouch Vita, contact InTouch Health and discontinue use of the system until the problem can be remedied.

The following tables contain the Manufacturer's declaration and additional information required by IEC 60601-1-2.

**WARNING:** Leakage current from interconnected electrical equipment may exceed safe levels. In order to maintain patient and user safety, it is important to interconnect only with devices in compliance with IEC 60601-1-1 requirements. It is the responsibility of the user to ensure that any interconnected equipment not supplied by InTouch Health maintains compliance with IEC 60601-1-1 requirements.

**Table 1: Electromagnetic Emissions**

<p>The InTouch Vita is intended for use in the electromagnetic environment specified below. The customer or the user of the InTouch Vita should assure that it is used in such an environment.</p>		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment</b>
RF Emissions CISPR 11	Group 1	The InTouch Vita uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The InTouch Vita is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations / Flicker Emissions IEC 61000-3-3	Complies	

## Table 2: Electromagnetic Immunity

The InTouch Vita system is intended for use in the electromagnetic environment specified below. The customer or the user of the InTouch Vita should assure that it is used in such an environment.

Immunity Test	EC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV Contact ±8 kV Air	±6 kV Contact ±8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient / Burst IEC 61000-4-4	±2 kV for Power Supply Lines ±1 kV for Input / Output Lines	±2 kV for Power Supply Lines ±1 kV for Input / Output Lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV Line(s) to Line(s) ±2 kV Line(s) to Earth	±1 kV Line(s) to Line(s) ±2 kV Line(s) to Earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Main power quality should be that of a typical commercial or hospital environment. If the user of the requires continued operation during power mains interruptions, it is recommended that the InTouch Vita be powered from an uninterruptible power supply (UPS) or a battery.

Immunity Test	EC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power frequency (50/ 60 Hz) Magnetic Field IEC 61000-4-8	3 A/ m	3 A/ m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
U <sub>T</sub> is the a.c. mains voltage prior to application of the test level.			

### Table 3: Electromagnetic Immunity

InTouch Vita is intended for use in the electromagnetic environment specified below. The customer or the user of an InTouch Vita should assure that it is used in such an environment.

<b>Immunity Test</b>	Conducted RF IEC 61000-4-6 Radiated FR IEC 61000-4-3
<b>EC 60601 Test Level</b>	3 Vrms 150 kHz to 80 MHz 3 V/ m 80 MHz to 2.5 GHz
<b>Compliance Level</b>	3 Vrms 150 kHz to 80 MHz 3 V/ m 80 MHz to 2.5 GHz
<b>Electromagnetic Environment - Guidance</b>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the InTouch Vita, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> <p><math>d = 1.2 \sqrt{P}</math>            80 MHz to 800 MHz            800 MHz to 2.5 GHz</p> <p><math>d = 1.2 \sqrt{P}</math>  <math>d = 2.3 \sqrt{P}</math></p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**NOTES:**

1. At 80 MHz and 800 MHz, the higher frequency range applies.
  2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.
- 
- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the InTouch Vita is used exceeds the applicable RF compliance level above, the InTouch Vita should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the InTouch Vita.
  - b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [3] V/m.

## Table 4: Recommended separation distances

### Recommended separation distances between portable and mobile RF communications equipment and the InTouch Vita.

The InTouch Vita is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the InTouch Vita can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the InTouch Vita as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

#### NOTES

1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Safety Warnings and Cautions

### WARNINGS

- InTouch Health requires all users to first be trained in the proper use of the InTouch Vita. The InTouch Vita stands approximately 5.5 ft. tall (168 cm) and weighs about 180 lb. (82 kg). An untrained operator could potentially bring about a collision, possibly causing damage or injury.
- The InTouch Vita is not MRI (Magnetic Resonance Imaging) safe, and therefore not MRI compatible. The InTouch Vita should not be used in locations where the presence of metal is controlled.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the InTouch Vita, including cables specified by InTouch Health. Otherwise, degradation of the performance of this equipment could result.
- Flammable Anesthetics: The InTouch Vita is not suitable for use in the presence of flammable anesthetic mixture with air, or in the presence of a flammable anesthetic mixtures with oxygen or nitrous oxide.
- The InTouch Vita must be plugged into a grounded “Hospital Grade” electrical outlet to minimize the risk of electrical shock during the battery charging cycle.
- Provide adequate ventilation. The InTouch Vita may overheat if powered on or plugged in and charging while stored for extended periods in an area without adequate ventilation.
- Leakage current from interconnected electrical equipment may exceed safe levels. In order to maintain patient and user safety, it is important to interconnect only with devices in compliance with IEC 60601-1-1 requirements. It is the responsibility of the user to ensure that any interconnected equipment not supplied by InTouch Health maintains compliance with IEC 60601-1-1 requirements.
- InTouch Patient Access Devices are designed to utilize the 802.11 ac/ a/ b/ g/ n communication protocol as well as the public Internet in order to achieve their intended purpose. Failures in either of these supporting systems could cause a complete loss of communication between the

Provider Access Software and the Patient Access Device. A tested backup method of communication should be available in the event network communication is lost.

- The video images transmitted to and displayed on the Patient Access Device and InTouch Provider Access Software may not contain all of the information in the original scene. Video information from the camera is captured, compressed, transmitted, and redisplayed re-motely at a different resolution. As a result, information in the original scene may be lost.
- Color reproduction in the transmitted video is not guaranteed. Color reproduction in a video system is a complicated combination of lighting, cameras, and display technology. It should not be assumed that the colors on the display are an exact replication of the actual colors in the scene.
- Avoid reliance upon video systems for diagnosis and video teleconferencing. Limitations including those outlined above mean video conferencing is not a perfect substitute for in-person interaction with a trained care provider.
- Clinical judgment and experience are required to review and interpret images and information transmitted via the Patient Access Device and Provider Access Software.
- Do not attempt to open or remove any parts of the InTouch Vita. To reduce the risk of electric shock, do not remove the cover. There are no user-serviceable components inside. Refer servicing and repair to qualified personnel only.

### CAUTIONS:

- InTouch Health does not support the addition of third party software to an InTouch Vita. Adding third party software (especially for video conferencing) to the computer can cause the InTouch Vita to malfunction. Please be advised to check with Technical Service PRIOR to installing any third party software.
- The InTouch Vita should be plugged in whenever it is possible so it is fully charged and ready for a consult.
- There are no user-serviceable components. Refer servicing and repair to qualified personnel only.
- DO NOT IMMERSE the InTouch Vita. DO NOT ALLOW any cleaning solution inside the InTouch Vita. Avoid excess solution which may enter the InTouch Vita through its openings.
- Keep the InTouch Vita free from moisture and extreme temperatures.
- InTouch Health has not performed safety and efficacy testing for many peripheral USB devices being used with the InTouch Vita. Customers must test and validate medical peripherals for their own use cases and environments
- Ensure external USB devices are disconnected prior to moving the InTouch Vita.

# Virtual Care System Overview

Health systems now view virtual care as an extension of their services, relying on a combination of software, hardware, networks, systems, and people to work together to deliver improved access, quality, and care to their patients.

Enabling healthcare's only integrated virtual care platform, InTouch Health powers virtual visits within your clinics, healthcare facilities, or even your patients' home for an integrated experience across a multitude of use cases, environments, and budgets.


Built on our cloud-based network, Solo™ by InTouch, our next generation software platform, is the backbone to delivering care anywhere, anytime. Solo by InTouch provides users with everything they need to streamline their telehealth needs for fast user adoption and program growth.

With Solo by InTouch as your virtual care foundation, connect people, healthcare systems, HIT systems, devices, and third-party software applications to enable telehealth solutions across the patient care journey. Whether launching a scheduled visit from your EMR, providing care for patients at home or in a clinic, Solo by InTouch provides a core set of tools to easily deliver virtual care and provide users with an immersive, patient-centric view of their clinical workflows.

## **Designed for healthcare, security, and reliability**

Our cloud-based, patented network ensures the industry's highest connection success rate and meets industry-defined requirements for protecting and securing sensitive healthcare information. Our downloadable and web-based platform allows users to access virtual care across a broad range of consumer and telehealth devices in any clinical setting.

Within a single platform, users are in control of their virtual experience with a configurable dashboard to meet the unique set of preferences for every user. Layer Solo by InTouch into your existing HIT system investments or use it as the foundation to close your workflow gaps with our software modules, creating a solution unique to you.



# InTouch Vita and Its Intended Application

The InTouch Vita is a clinical telecommunications tool that provides a means of transmitting and receiving data. InTouch Vita may also be used in conjunction with 510 (k)-cleared devices that transmit patient biometric data including vital signs information. The InTouch Vita transmits and receives information over an internet connection between patients and health professionals. The InTouch Vita can be used for telecommunications in both low and high acuity clinical use cases environments, e.g. pre-, peri-operative and post-surgical, cardiovascular, neurological, prenatal, psychological and critical care assessments and examinations.

## InTouch Vita Indications for Use

The InTouch Vita is used in telecommunications for a variety of clinical applications. These applications include, but are not limited to: Intensive Care (ICU), cardiovascular, neurological, behavioral health, neonatal, etc., use cases. The InTouch Vita does not diagnose health care conditions but instead provides the means for the remote physician to administer his or her assessments and examinations. Clinical judgment and experience are required to review and interpret the information transmitted.

# HIPAA Policies and Procedures

As a business associate, InTouch Health is subject to the compliance of the law under §164.308 (Administrative Safeguards), under §164.310 (Physical Safeguards), and under §164.312 (Technical Safeguards) in that we maintain and transmit health information in electronic form in connection with transactions performed by the customer (covered entity). The policy of this organization is to ensure, to the greatest extent possible, that Protected Health Information (PHI) is not intentionally or unintentionally used or disclosed in a manner that would violate the HIPAA Privacy Rule or any other federal or state regulations governing confidentiality and privacy of health information.

There are a number of safeguards implemented into the telehealth system that ensure that it complies with the latest HIPAA regulations. One of the key requirements is InTouch Health's ongoing implementation and updating of its HIPAA security policies and procedures to ensure for the availability, security, and privacy of telehealth connections and ePHI (electronic protected health information). InTouch Health maintains a policy to ensure workforce HIPAA compliance and training. InTouch Health additionally maintains HIPAA security policies and procedures, a data destruction policy, and security incident response procedures.

## Recommended Guidelines for HIPAA Compliance

The telehealth system allows hospitals and medical professionals to remain in compliance with HIPAA privacy regulations. InTouch Health is providing the following information to assist with HIPAA privacy regulations as they pertain to the telehealth system.

HIPAA requires all healthcare organizations to have policies and procedures, and the guidelines below may not cover all situations as they apply to a specific organization. Further, from time to time, automatic software upgrades may be downloaded which may contain new features. InTouch Health will inform users of significant features as they are added, and the impact of these features and how they may affect HIPAA policies, procedures, and safeguards should be considered.

## Access to Provider Access Software

The computer using the Provider Access Software should be placed in a location that is only accessible to individuals who have authorized access to Protected Health Information (PHI). It is recommended that Provider Access Software access is password protected via a Windows or iOS user account.

Only authorized users should have passwords, and users should safeguard passwords according to hospital policies and procedures. Passwords should be treated as highly confidential information. If you believe your password may have been compromised, it should be changed as soon as possible. Change your password by clicking on the "Forgot Password" link on the login screen of the InTouch Provider Access Software.

The Auto Logout feature is set to log out of the InTouch Provider Access Software when the system is inactive for 30 minutes. Also, all users should be trained to log out of Windows, iOS or the Virtual Private Network (VPN), when away from the system for any period of time. This is important for security reasons, so that any person attempting access to the Provider Access Software will be required to enter a password for secure access.

## Discussion and Display of PHI

From time to time a physician will likely engage in remote communications with patients and medical staff in which patient information (records, images and video) will be discussed or displayed. In general, the same care should be exercised as though the physician were physically present. For example:

- Use Head rotation to look around and see who else is nearby and might see or hear the sensitive information, and use appropriate discretion.
- Use the microphone mute button when conversing with someone alongside the InTouch Provider Access Software to avoid the inadvertent conferencing of patient-related conversation.
- The InTouch Provider Access Software screen should be positioned to point away from public areas, so as not to be visible to passersby.

## Stored Images and Video

By default when saved, all captured images and video files are stored as encrypted files viewable only by the Provider Access Software user who captured them. All files whether encrypted or not should also be saved in the user's InTouch Health Media Vault (folder) to provide added protection.

For convenience these files may be saved in common formats, e.g., JPEG for still images. These files are no longer encrypted and therefore are viewable by any user who can access them. As such, there are a few recommended techniques for safeguarding PHI contained in these images and video:

- Ensure all personnel who have access to the Provider Access Software also have full permission to access stored images and videos under the hospital's policies and

procedures;

- Make sure to store captured images and videos only on removable media (e.g., recordable CD-ROMs) which can be taken with each user or on secure network drives;
- Do not save to disk any captured images and video clips. Use these images and video segments only while logged in for a particular session.

## Disclosure of PHI

If the physician plans to transmit or copy stored images or video to other individuals or organizations, e.g., to a healthcare operator, the physician needs to abide by standard HIPAA codes governing who may receive PHI and under what conditions. The hospital's HIPAA compliance officer should be consulted for details.

# InTouch Vita Overview



Chest Display	Touch-screen used for local control of Vita.
Collision Avoidance Sensors	InTouch Vita's collision avoidance sensors consist of 3D near-IR depth sensors and sonars to provide for detection of obstructions.
Directional Microphone	Captures audio for playback on the Provider Access Software.
Docking Charger Interface	Used to connect to Dock (allows charging and automated Docking).
Emergency Stop Button	Stops Vita's (base) movements by disengaging its motors. Also, used to disengage wheels, when needing to quickly move Vita manually. See "Emergency Stop Button".
Head Speakers	Plays audio from Provider Access microphone on Vita.
Laser Pointer	FDA Class II Laser Pointer fixed to the Vita pan/tilt Head. Avoid exposure. Laser radiation is emitted from this aperture.
LED Indicator Lights	Vita status indicators. Refer to "Vita Status Indicator Lights".
Power Button/ Power Cord (inside cord bay door)	The Vita power button turns on the computer and display. Power cord used to recharge Vita without the Dock.
Remote Display	Displays remote user's face (video from the Provider Access) or recorded pictures and video from the Provider Access on Vita.
Rocker Switch (Main On/ Off Switch)	Removes power to the VITA computer systems.
Storage/ Expansion Bays	An expansion bay is located on each side of the Vita. One side contains the Stethoscope and Headset, the other side contains the Privacy Handset and the Input/ Output ports for connecting approved devices to Vita.

Virtually There Cameras	Capture remote video for viewing at the Provider Access.
Volume Control Knobs	Used to adjust volume. See “Volume Adjustment Knobs”.

The Vita utilizes a Class II Laser Pointer which complies with 21 CFR Chapter 1, sub-chapter j. Maximum Laser radiation output is less than one milliwatt at a wavelength of 635 nanometers.



**CAUTION:** Dazzle, flash-blindness, and after-images may be caused by a beam from a Class II laser product, particularly under low ambient light conditions. This may have indirect general safety implications resulting from temporary disturbance of vision or from startle reactions. Such visual disturbances could be of particular concern connected with performing safety-critical operations. Users should not stare at the beam and perform active protective reactions by moving the head or closing the eyes to avoid continued intrabeam viewing.

**WARNINGS**

- The Vita is not MRI safe nor MRI compatible and should only be used in locations where the presence of metal is not controlled.
- Flammable Anesthetics: The Vita is not suitable for use in the presence of flammable anesthetic mixture with air, or in the presence of a flammable anesthetic mixture with oxygen or nitrous oxide.

# InTouch Vita Basics

## Charging the InTouch Vita

- InTouch Vita will charge to 80% in two hours.
- InTouch Vita will charge to 100% in six hours.

**NOTE:** Keep Vita docked whenever it is not in use.

The InTouch Vita can be sent to its Dock from the Provider Access Software or by hospital personnel using the Touch-Screen Chest Display. See AutoDRIVE.

If Dock is not available for any reason and InTouch Vita needs to be charged, plug the Vita into a standard grounded outlet.

**NOTE:** The power plug has an internal light which lights when plugged into a live outlet. Ensure it is lit when charging with the power plug.

**CAUTION:** The InTouch Vita contains high capacity rechargeable, Lithium-Ion batteries. The InTouch Vita should be plugged in when not in use so it is fully charged and ready for the next consult and to avoid deep discharge cycles that can shorten the battery's useful life. Other than keeping the batteries charged by keeping the Vita plugged in, no user maintenance of the batteries is required.

## Charging Status Indicator Lights

The InTouch Vita Dock has an LED indicator light located on the top of the Dock. In addition, the power plug contains a light which turn on when plugged into a live AC outlet.

Dock	LED State
No AC power applied	Off
AC power applied, not mated to robot	On
AC power applied, mated to robot	On - blinking
InTouch Vita Power Plug	Light State
Plugged in and powered	On
Unplugged or not powered	Off



## Vita Power Controls

The InTouch Vita has two power controls:

- the Rocker Switch (main On/ Off switch), that removes power to the InTouch Vita computer systems; and
- the InTouch Vita Computer Power Button.

The Rocker Switch is located just below the skirt on the base, above the Docking Charger Interface.

The Rocker Switch, which remains in either the Off or On position, controls the power to the InTouch Vita. If the InTouch Vita is neither plugged in nor docked, the power is provided from the battery. The Vita Rocker Switch should stay On, except during storage, service, transport, or to avoid depleting the battery.

### Power On Sequence

1. Switch the Rocker Switch to the On position.
2. Plug the InTouch Vita in at its docking station or plug the power cord into an approved wall socket.
3. Turn On the power to the computer:

Press and hold in the Computer Power Button for about half of a second.


Check to ensure the display powers up successfully.

After a period of self-testing (1-2 minutes), the screen saver should appear on the display.

**NOTE:** Keep the InTouch Vita docked (or plugged in) and left powered on at all times.

#### CAUTIONS:

- To ensure system readiness, connectivity, and charged batteries, power on the InTouch Vita at least two hours before its intended use.
- The InTouch Vita checks for and installs any available and scheduled updates automatically when it is docked.



## Power Off Sequence


**CAUTION:** When powering down the Vita for any reason, always ensure that the Computer Power Button is powered off first and that the Vita Display screen turns off (turns black) before turning off the Rocker Switch. It will take approximately 30-45 seconds for the Display screen to turn off.

1. Turn Off the Power to the computer:
  - Press the Computer Power Button.
  - Check to ensure that the display screen turns off (turns black).

### NOTES:

- It will take approximately 30-45 seconds for the screen to turn off.
- If the computer does not shut down following the above step, you can force a shut down by holding the Computer Power Button in for several seconds until the Vita Display screen turns black.

2. Turn Off the Rocker Switch.



## Vita Emergency Stop

### To disable Vita movement:

Press the red Emergency Stop Button.

- The Robot computer will stay On, but the Vita cannot be driven as its base motors are disengaged.
- The Robot can now be pushed where needed.

### To enable Vita movement again:

- Twist the red Emergency Stop Button clockwise approximately 1/4 turn until it pops out.



## Vita Volume Adjustment Knobs

The Vita has two volume controls: one on each side of the Touch-screen Display.

- Turn either dial to control the volume of the Vita speakers.




## Vita Head and Chest Display



- InTouch Vita can be accessed by a remote provider using InTouch Health's software. A provider can connect to the Vita without interaction on the Vita, as a powered-on Vita automatically accepts a connection request from any authorized provider who knows the device's serial number or Care Location's name (available on Diagnostic Information screen).
- When Vita is actively connected to Provider Access Software, the video from the Provider Access Software will be displayed on the Vita Head Display and the name of the person in the session will be shown on the Chest Display.
- When the Vita is NOT actively connected to Provider Access Software, it is considered Idle. When Vita is idle, an animation will be displayed. Screen savers selected by the customer may optionally be shown; to do so, consult your InTouch Health representative or Technical Support.

## Vita Touch-Screen Chest Display


The Vita Touch-Screen Chest Display provides In Session information and local controls of the Vita. Whenever the Vita is powered on, you can tap the Chest Display to interact.



Exit	Puts Vita into Idle state, with animation on the screens.
AutoDRIVE	Choose a pre-defined destination for the Vita to travel to.
Camera	Camera Preview mode allows the local user to preview the main camera image, and that of any attached video peripheral.
Settings	Access device settings

Info	Diagnostic Information, such as robot serial number, battery charge, Wi-Fi signal strength, Robot software version, InTouch Health's Technical Support contact info, and other data.
------	--

## Vita AutoDRIVE Screen



To send Vita to a predefined location, tap the AutoDRIVE button, then tap the desired destination location button.

While in motion, tap the Chest Display to halt Vita.

- When stopped, the Chest Display shows Resume and Cancel buttons.
- Tap the Resume button to allow Vita to continue on to its destination. If not resumed or canceled, the navigation will time out.

The Vita's Chest Display will report when it reaches its destination or if it fails to reach its destination.


## Docking the Vita

To initiate automated docking, when it is available, tap the AutoDRIVE button, then the Dock button

**NOTE:** Dock will always appear at the top of the list.

## Change AutoDRIVE Map

The Vita has Multi- Map capabilities that allow it to utilize AutoDRIVE on different floors of your hospital. During installation, please ask the Product Specialist to map each floor in which the remote physician will use AutoDRIVE.



Follow these steps to switch floor maps on the Vita:

**CAUTION:** Make sure to check destinations before and after the map change.

1. Tap anywhere on the chest display of the Vita. Tap Settings on the Tool bar.
2. Tap on Change Map.
3. Tap on the new floor under Available Maps that corresponds to the new location of the Vita.
4. The Vita will prompt you to confirm that you want to switch maps. Click Yes to confirm.
5. The Vita will display a message alert that reminds you to dock the Vita before use.
6. The map should have successfully been changed. The highlighted floor name under Available Maps is the current map on the Vita.

**NOTE:** Map can also be switched during a live session using steps 1-6 above. There will be a confirmation prompt, and Message Mode will not appear. The remote provider will not have AutoDRIVE capability if the Vita is not properly docked.


Set BLG2-Vita Map-QA-Patrol as the current map?

Set

Cancel


## Vita Settings Screen

Tap  on the InTouch Vita's chest display to show the available settings.




## InTouch Vita Wi-Fi Connection Setup

1. Select **Wifi Setup**  from the **Settings** screen .
2. Select the desired secure, clinical wireless SSID network.



**CAUTION:** Guest, Staff, and VIP type networks are not suitable for connection of a medical device.

3. Enter the network security key.



4. Use the **Information Icon** (  ) on the lower right corner of the screen and ensure that the status is **Ready** before attempting to make a connection.

**NOTES:**

- The Diagnostic Information page will display current WiFi status.
- InTouch Vita will automatically connect to known Wi-Fi profiles. For stored WiFi profiles, the password can be changed by selecting it and editing the password.

**Connectivity Details**

The InTouch Vita can connect to the Internet using cellular broadband through MiFi or other devices capable of creating WiFi hotspots. The cellular service provider and/or signal strength may change from area to area. Customers should contact InTouch Health Technical Service for further information.

- Wireless Network Signal Strength Indicator is located on the Diagnostic Information page
  - If Indicator states Fair user may experience video problems.
  - When Indicator states Weak the Robot may lose the connection—it is best to return to a stronger signal area immediately.

Ready 

### Diagnostic Information

#### InTouch Vita Ready

- Battery Charge ↑ 99%
- Bus Voltage 36.09
- Serial Number Robot-9369
- Location Demo-Vita 2
- Software Version 10.39.89
- Build Date August 06 2019
- Autodrive Status Available
- Base Status Connected
- Stethoscope Type PCP

#### Network Details

- Connection Type WiFi
- WiFi SSID TIC
- WiFi Strength ↑ 100%
- External IP 64.29.227.117
- Internal IP 192.168.12.17

#### Server Connectivity Status

- SIP ● Registered
- FMS ● Registered
- WebRTC ● Registered


 Customer and Technical Services  
(877) 484-9119

## Device Cleaning Mode

The InTouch Vita has a screen-cleaning mode, so that no functionality is inadvertently called upon when cleaning.

- See ["Cleaning Instructions" on page 50](#) for more details.


Tap Settings  > Cleaning Mode , to disable the screen for 30-seconds to allow cleaning.



## Date and Time Settings


Tap Settings  then Date and Time  to adjust how the date and time are displayed on the InTouch Vita during a session.

Tap  to return to the Settings screen.



## Device Settings - More Options


From the **Settings** screen, select the  to display additional device settings.



## Vita Change Map



1. Tap the  from the Settings screen to select different, preprogrammed maps.
2. Scroll through the available maps to locate the appropriate one.
3. Tap the appropriate map to select it.
4. Tap  to return to the idle screen.




## Vita Edit Badge


This feature allows the information found on the Chest Display name badge.




1. Tap  from the Settings screen to edit the information on the InTouch Vita's Chest Display's badge.



2. Tap the field to change the text.
3. Use the onscreen keyboard to enter the revised text.



4. Tap  to accept the new text.
5. Select a different icon, if available, by tapping the left or right arrows.
6. Tap  to return to the home screen.

## Vita Status Indicator Lights

InTouch Vita has LED light strips that indicate its operational status and condition. Lights are located on all four sides under the base and on the left and right sides just above the base. The Head and Chest Displays will also display the color. The meaning of the indicator lights is described in the following table. See “Vita Anatomy and Components” for the location of the LED indicator lights.

Vita Status	LED Color, Animation	Example
Start Up	<ul style="list-style-type: none"> <li>• Grey Base lights</li> <li>• No Head or Chest display</li> </ul>	
Idle (not charging)	<ul style="list-style-type: none"> <li>• White/ Gray color pallet</li> <li>• IDLE badge displayed</li> <li>• Badge displays greeting</li> </ul>	
Idle (charging)	<ul style="list-style-type: none"> <li>• White/ Gray color pallet</li> <li>• Idle badge displayed</li> <li>• Badge displays greeting - Base lights pulse</li> </ul>	
Message—See “Message Mode”	<ul style="list-style-type: none"> <li>• Magenta color pallet</li> <li>• Message badge displayed</li> <li>• Notification shows error</li> </ul>	

Vita Status	LED Color, Animation	Example
Mission—AutoDRIVE command from the Chest Display other than Dock.	<ul style="list-style-type: none"> <li>• Blue color pallet</li> <li>• Mission badge displayed</li> <li>• Badge displays destination</li> </ul>	
Session	<ul style="list-style-type: none"> <li>• ITH Teal color pallet</li> <li>• Session badge displayed</li> <li>• Badge displays User name</li> </ul>	
Fault	<ul style="list-style-type: none"> <li>• Yellow/ previous color pallet</li> <li>• Previous badge displayed</li> <li>• Notification shows error - Yellow LED light</li> </ul>	

**NOTE:** If InTouch Vita is going to its Dock and is not currently in a session, the LED colors will be static white. This can occur as a result of InTouch Vita being sent to the Dock from the Chest Display, a disconnected session, or losing its Wi-Fi signal, or because of an idle or low battery timeout.

# Cleaning and Maintenance

## Cleaning Instructions

Clean the outer surfaces of the Device when visibly soiled or after contact with any contaminants. All surfaces, such as display monitors or sensor windows may be disinfected using the following procedure. Use a commercial LCD screen cleaner to prevent craze, staining, or discoloration of the display monitors and use optical lens cleaners to clean the camera lenses.

**NOTE:** For devices with touch screens, go to Settings > Cleaning Mode, to temporarily disable the touch interface for 30-seconds to enable cleaning. See ["Device Cleaning Mode" on page 42](#)

1. Power off and unplug the device.
2. Soak a lint-free cloth in a hospital grade disinfectant solution of sodium hypochlorite 6.15%, e.g., dilution 1:500 (1/4 oz. per gallon water) and wring out the cloth so that drips do not appear when wiping surfaces.
3. Wipe surfaces that have become soiled or contaminated. Avoid applying excess solution which may enter the device through its openings.
4. Allow to air dry.

## Approved Disinfectants

The following disinfectants have been tested for compatibility with InTouch Health devices:

**NOTE:** the following agents are for the exterior, not the lens or screen.

- OxyCide
- Ethyl Alcohol
- Isopropol Alcohol
- Sodium Hypochlorite (5.25%- 6.15% household bleach diluted 1:500 provides >100 ppm available chlorine).
- Iodophor Germicidal Detergent Solution (follow product label for use dilution).
- Quaternary Ammonium Germicidal Detergent Solution (follow product label for use dilution).

### **WARNINGS:**

- Do not attempt to open or remove any parts of the InTouch Vita.
- Do not remove any covers to reduce the risk of electric shock. There are no user-serviceable components inside.
- Refer servicing and repair to qualified personnel only.
- Wear safety glasses when handling solution prior to dilution.
- Wear rubber or nitrile gloves, if in contact with liquid.
- Avoid contact with eyes, skin and clothing.
- Wash hands after cleaning device.
- Do not wear product-contaminated clothing for prolonged periods.
- Always follow manufacturer's instructions on product labels when mixing chemicals.

### CAUTIONS:

- **DO NOT USE** phenolic germicidal detergent solutions on any parts of the device. Contact InTouch Health Technical Support for approved cleaning solutions.
- Severe contamination, especially of the undercarriage may require some disassembly and this should only be done by an InTouch Health representative.
- **DO NOT IMMERSE** the device.
- **DO NOT ALLOW** any cleaning solution inside the device.
- Keep the device from moisture and extreme temperatures.
- Avoid excess solution that could enter the device through its openings.

## Cleaning Sensor Windows (Vita)

Keep the sensor windows clean to ensure dirt or foreign particles are not interpreted as objects to avoid when the device is moving.

- Clean all three windows of the waist sensor.
- Clean the entire 270° extent of the sensor located under the skirt at the front of the device.



## Maintenance and Inspection

### InTouch Vita Maintenance

The InTouch Vita contains no user serviceable parts and requires no maintenance. For further information regarding preventive maintenance, maintenance or assistance with troubleshooting, customers should contact InTouch Health Technical Service.


### Vita Inspecting the ESD Cable

The Electric Static Discharge (ESD) cable is located on the underside of the Vita to dissipate any static charges that arise. The cable should be checked periodically to ensure that it is present and remains in contact with the floor.

**CAUTION:** Do not tilt Vita onto its bumpers or attempt to place Vita on its back. The bumpers located on the bottom edge of Vita may be damaged.

- Tilt the Vita toward the back on its roller wheels so that the static cable can be viewed.
- As shown in the photograph, the cable should be relatively straight and not kinked.
- Verify that the cable reaches the floor when the Vita is upright.

Contact InTouch Health at **(877) 484-9119** if the ESD cable appears damaged or is missing.



# Vita Stethoscope

## Vita Stethoscope Components

- Stethoscope assembly (pre-mounted on the Vita)
- Stethoscope chest piece
- Vita headset

## Vita Stethoscope Operation

On the Vita side, the Stethoscope chest piece is applied to a patient by hospital staff following the physician's directions (received through normal Vita audio).

**CAUTION:** Do not pull on the cables, and do not pull the Vita using the cables as cables may break, especially near the junction points with the Vita.

1. Ensure the Vita is within 6 feet of the patient exam site.
2. Put on the Stethoscope headset (if desired).
3. Put on gloves.
4. Remove the Stethoscope chest piece and disinfect it with a disposable sanitary wipe.
5. Expose the patient area and apply the chest piece as directed by the physician.
6. Disinfect the chest piece and return it to its bracket.
7. Remove gloves.
8. Remove the headset and return it to its bracket.
9. Close the Storage bay door.



## Vita Stethoscope Headset

Using the Vita headset is optional, but it may be used to assist with placement of the chest piece on the patient or as an aid for teaching and mentoring. The Vita headset allows the nurse, resident, fellow, etc. to listen to the Stethoscope while still being able to hear the normal Vita audio and other sounds in the environment.

## Technical Service—Stethoscope

For any issues with the PCP Stethoscope on the Vita, contact InTouch Health Technical Service. Questions about the PCP Stethoscope controls can be answered by reference

to the PCP Stethoscope manual shipped with the Stethoscope kits. If a chest piece or headset cable breaks, contact InTouch Health Technical Service for assistance.

**CAUTION:** Do not pull on the cables, and do not pull the Vita using the cables as cables may break, especially near the junction points with the Vita.

# Privacy Handset

The Privacy Handset is for use in conversations that require privacy or understandability in a loud environment.

When the Privacy Handset is enabled from the Provider Access Software, the standard Patient Access Device audio system is disabled. This means the Patient Access Device microphone is not active, and the Patient Access Device speaker is not active.

When the Privacy Handset is disabled, the audio system is transferred from the phone, back to the standard Patient Access Device audio system.

It is important to return the Privacy Handset to the hook on the Patient Access Device. Ask someone on the Patient Access Device-side to return the Privacy Handset to the hook if necessary.



**NOTE:** Always place the Privacy Handset on the hook after use.

# InTouch Vita Aux Video Ports

The InTouch Vita's expansion bay includes both video input ports and a VGA output port. The USB ports can also serve to capture video from approved USB cameras. This allows the remote physician to view multiple video inputs through the Provider Access Software.

The VGA output allows the video from the Provider Access (shown on the InTouch Vita remote display) to be shown on a different screen or on a projector.

## WARNINGS:

- The video images transmitted to and displayed on the InTouch Vita and Provider Access may not contain all of the information in the original scene. Video information from the camera is captured, compressed, transmitted, and redisplayed remotely at a different resolution. As a result information in the original scene may be lost.
- Color reproduction in the transmitted video is not guaranteed. It should not be assumed that the colors on the display are an exact replication of the actual colors in the scene.
- Clinical judgment and experience are required to review and interpret images and information transmitted via the InTouch Vita and Provider Access.

**CAUTION:** Always push the Emergency Stop button in when connecting video equipment to the Auxiliary Video Input. Disconnect the equipment when ending a session and before pulling the Emergency Stop button out to enable the InTouch Vita's (base) movements. The InTouch Vita has Autonomous drive modes that will be commanded when it is not in session. For example, the InTouch Vita will return to its Dock when its battery charge is running low. The InTouch Vita and/or attached equipment may be damaged if the InTouch Vita is driven with equipment connected to its Inputs.



# InTouch Vita Physical Specifications



## Vita Technical Specifications

<b>Performance</b>		
Head:	Pan range: +/- 170° max Pan speed: 90°/sec max	Tilt range: 127° max Tilt speed: 90°/sec max
Base:	Maximum speed: 3.36 mph (forward) Locomotion: holonomic drive system	
Sensors:	Multiple center and base-mounted cameras and sensors	FDA Class 1 Laser sensors detect obstructions up to 98 ft. (30m) within a 270 degree range centered about the front of the robot.
Audio:	Microphone: Mono, directional (hyper-cardioid), 50Hz-19kHz, Speakers: two 5 W mono speakers in the head and a sub-woofer in the base. Capable of 100dB sound (1KHz @1m).	16 kHz sampling rate, 16-bit audio
Video:	Camera: 120X equivalent zoom, remote zoom & focus Video: 30 fps, 648x480 px resolution, 24-bit color Display: 15 in LCD, 1024x768 Touch-screen: 8.4 in. LCD, 1024x760	Display: 15" LCD, 1024x768 px, 400 Nits
Wireless Network:	802.11 a, b, g, or n	
Battery Life:	4-5 hrs (depending on usage)	
Charging Time:	2 hours from 0% to 80% (6 hours to 100%)	
Driving Range: (distance)	Related to the range of the wireless network in use	

## Vita System Input Power Requirements

	North America	International
Voltage	115 VAC	230 VAC
Frequency	60 Hz	50 Hz - 60 Hz
Current	10 Amps	5 Amps

## Vita Classification

- Internally Powered / Class I while in Charging Mode
- Type B applied parts - Headset and **Handset**
- Type BF applied parts - Stethoscope

## Vita Environmental Specifications

### Operating

Designed to operate in an indoor environment suitable for human personnel.  
(10° to 30° C, 30 to 75% RH, 700 hPa to 1,060 hPa)

### Non-operating/Storage/Transport

Designed to travel to installations in commercial and cargo airliners and standard ground transportation. (-40° to +70° C, 10 to 100% RH, 500 hPa to 1,060 hPa)

# Network Configuration

**NOTE:** For the best performance from the device in terms of Device Optimization and connection success, utilize the information found in these two documents:

- MB-15513 - Network Configuration for InTouch Devices
- MB-14011 - InTouch Telehealth Network Specifications

Contact your InTouch Health representative for copies of these documents.

The InTouch Telehealth System is comprised of a Remote Presence Patient Access Device and a minimum of one Provider Access Software Device. The Provider Access and Patient Access Devices are linked via the Internet over a secure connection.

The optimal configuration is determined during installation.

## Configuring the InTouch Patient Access Device Wireless Connection

The InTouch Vita Control Core uses a Windows computing environment and a wireless network card.

For the InTouch Vita, basic Wireless Network connections can be made using the touch-screen, under Settings.

## Non-Overlapping Channels

In order to achieve a smooth transition from one wireless access point (WAP) to the next it's important to configure each WAP on a non-overlapping channel.

## Transmitting Power

The WAPs may be configured to provide the wireless signal at different power transmission levels. Setting the WAP to the maximum power transmission will deliver the maximum coverage area.

## Interference


If WAPs are co-located in the same environment, radio frequency interference may be generated. Too many WAPs transmitting on overlapping channels may also degrade the wireless signal quality.

WAPs placed too close to one another may also produce RF congestion. In this case, the WAP transmission power should be reduced; therefore, reducing the coverage area and limiting the overlap between adjacent WAPs.

## Security Options

Each wireless network must be configured with security to prevent unauthorized access to the network. The ADU provides multiple features to configure the InTouch Vita to access as well as secure the wireless network. Domain membership is not supported, but all other current security configurations can be configured. WPA2/ AES- PSK is preferred.

# Network Installation



# Network Requirements

## Provider Access User Authentication

Users launch the Provider Access Software from their desktops and login to begin consults. Users have a unique username and password for their Provider Access accounts, which is created at registration.

Enterprise login allows users to use their company credentials to login to the Provider Access Software. If your system is configured to use Enterprise login, your system may automatically login to the Provider Access Software.

Users should make sure to comply to HIPAA standards by ensuring that they log off when away from their computers for extended periods of time. Users must log into Provider Access Software using their InTouch username and password. After a period of inactivity, users are logged out automatically as a security feature.

## Firewall Requirements

The InTouch Telehealth System uses bidirectional communication under TCP and/or UDP. The InTouch Telehealth Network consists of connections made through either our SharedComm or SIP servers. For optimal connections using SharedComm, Patient Access Devices and Provider Access Software require outgoing UDP access on ports 9000-9101 with reflexive UDP access ('UDP Replies') enabled. (Note: UDP replies are enabled by default on most firewalls). For optimal connections using SIP, Patient Access Devices and Provider Access Software require outgoing UDP access on ports 35000-35500 with reflexive UDP access ('UDP Replies') enabled.

For training and support purposes (including software upgrades), InTouch Health routinely makes use of remote desktop applications (Kaseya, GoToAssist & GoToMyPC). InTouch Health requires access to all Patient Access Devices and Provider Access Software via one of these applications.

For a detailed list of IP addresses and ports to white list, please refer to the document: Network Configuration for InTouch Devices (MB-15513). Additional modifications may be necessary for use with a Web Filter and/ or Stateful Packet Inspection.

## Video Information

**Frame rate:** Video is captured at 30 frames per second but can be reduced for low bandwidth connections.

**Codec (video and audio compression):** InTouch Health uses the standards-based H.264 AVC codec for video and the Opus or Speex codec for audio.

Dynamic Video Quality provides the ability to dynamically adjust resolution and video quality during a live session without user interaction. Advanced users can specify preferences for adjusting resolution. Video quality depends on factors such as robot motion, available bandwidth, and user preference.

## Bandwidth Requirements

The ideal bandwidth required is 700 kbps in both directions from any Provider Access Software or Patient Access Device. For Provider Access Software located in homes, lower bandwidths such as home cable broadband can be configured with good performance effectively utilizing 300 kbps. For installations where higher audio and video quality is desired, higher bandwidths above 700 kbps can be allocated.

For HD video, your Provider Access Device must be configured to allow 2000-3000 kbps.

## Line Quality Requirements

Network performance is critical to maintaining a responsive Provider Access Software to Patient Access Device session. Metrics cover a range of network characteristics which impact delivery of complete correct data in the proper order in a timely fashion. InTouch Health runs tests using proprietary and third party software tools to determine if a broadband connection meets a sufficient level of network performance to maintain a session. InTouch Health can provide these tools to customers upon request.

**NOTE:** During any particular session, quality may be degraded or the session may be disconnected if the network performance limits described below are exceeded, even though performance measurements were within limits at another point in time.

There are five important network characteristics affecting connectivity:

**Data rate:** A connection must have the required up-stream and downstream bandwidth, as discussed above.

**Latency (delay):** Average network latency on a connection should not exceed 300ms.

**Maximum Transfer Unit:** The Maximum Transfer Unit (MTU) must not be set below 1400 bytes.

**Reliability:** A connection must be reliable, without significant packet loss. A connection should experience no more than 3% packet loss.

**Jitter:** Jitter is variability in latency. Jitter on a connection should not exceed +/- 50ms during 95% of the duration of a session.

**NOTE:** InTouch Health's software is fully capable of dealing with the normal variability of data over the Internet. It is the quality of the endpoint connections which is critical and must be tested.

## Wireless Network Requirements

The InTouch Health System is compatible with 802.11 ac, a, b, g, and n protocols. The Maximum Handoff Threshold time must be less than 150ms.

In environments which experience network congestion, the InTouch Health application requires Quality of Service (QoS) or priority of traffic to ensure a successful connection.

## **Satellite Networks**

The network characteristics detailed above (bandwidth requirements, packet loss, jitter, and MTU) are strongly recommended to achieve an audio/ video session of functional quality over a satellite network. The one notable exception is the expected latency typical of satellite networks.

The InTouch Health System can maintain an audio/ video session of functional quality with latency up to 900 ms if all other network characteristics are met. Please note that this delay will be evident on both sides of any audio/ video/ command communication as is typical of satellite networks.

If utilizing a satellite network with latency above 600 ms, the delay in drive commands may hinder the operator from maintaining safe control over the movement of mobile devices. InTouch Health therefore does not recommend utilizing mobile devices, such as the InTouch 7i, on a satellite network.

## **Encryption**

The InTouch Health System incorporates encryption methodology utilizing a combination of RSA public/ private key and 256-bit AES symmetric cryptography. The following is a brief summary:

Each time a Provider Access session is initiated, a symmetric key is created using AES 256-bit cryptography. The encrypted data is then transmitted using RSA 4096-bit public-private key cryptography. PHI and other sensitive health information is further secured using SSL/ TLS and other different forms of authentication.

## **Virus Protection**

TrendMicro's OfficeScan is installed on every system. This software automatically updates as soon as new virus definitions are available. InTouch Health staff monitor software updates as they become available. InTouch Health staff install all necessary security updates on Patient Access Devices.

# Contact Information

24 /7 Live Technical Support

(877) 4 84 -9119

24 /7 Live Remote Technical Support & Live Chat

<https://intouchhealth.com/contact-us/>

E-mail Support

[ITHSupport@intouchhealth.com](mailto:ITHSupport@intouchhealth.com)

Website

[www.InTouchHealth.com](http://www.InTouchHealth.com)

InTouch Health User Manuals

<https://intouchhealth.com/manuals/>

Sales & Product Demos

(805) 562-8686



 InTouch Technologies, Inc.

7402 Hollister Avenue, Goleta, CA 93117

Ph: 805.562.8686 • Fax: 805.562.8663

[www.intouchhealth.com](http://www.intouchhealth.com)

